

Lichtstabilisatoren		
Parameter	CAS	Bestimmungsgrenze in µg/L
Tinuvin 234	70321-86-7	100
Tinuvin 312	23949-66-8	10
Tinuvin 326	3896-11-5	100
Tinuvin 622	65447-77-0	100
Tinuvin 770	52829-07-9	300
Uvasorb HA 88	136504-96-6	250
Uvinul 3030	178671-58-4	10
Uvinul 4050H	124172-53-8	2
Chimasorb 2020	192268-64-7	100
Chimasorb 81	1843-05-6	10
Chimassorb 2020	192268-64-7	100
Chimassorb 944	71878-19-8	100
Hostavin N30	202483-55-4	100
N,N-Bis(2,2,6,6-tetramethylpiperidin-4-yl) isophthalamide	42774-15-2	10